

System Testing

Sugar Analysis

CSE 115A Sugar Analysis Team:

Nooran Salim, Nikita Thumma, Pranav Nampoothiri, Justin Morales, Nicole Ng

Sprint One:

A. User Story One: As a product user, I need to be notified when I should test my blood sugar levels in order to maintain my health.

a. Did not finish in this sprint.

B. User Story Two: As a product user, I need the notification to transfer me to the homepage of the app so I can view my blood sugar check schedule.

- **Scenario: (test done by Nooran Salim)**

1. Start Sugar Analysis; opens to personal information page
 - a. Select Name
 - i. Textbox to <Enter Name> (Example: John)
 - b. Select Height
 - i. Textbox to <Enter Height> (Example: 68)
 - c. Select Weight
 - i. Textbox to <Enter Weight> (Example: 150)
 - d. Select Age
 - i. Textbox to <Enter Age> (Example: 21)
2. Home page has list view with time to input blood sugar
 - a. Before Breakfast
 - b. After Breakfast
 - c. Before Lunch
 - d. After Lunch
 - e. Before Dinner
 - f. After Dinner
 - g. Before Workout
 - h. After Workout

Sprint Two:

C. User Story One: As a product user, I need to check my logs stored in a database so I can monitor my progress.

a. Did not finish in this sprint

D. User Story Two: As a product user, I need the notification to transfer me to the homepage of the app so I can view my blood sugar check schedule.

a. Did not finish in this sprint

E. User Story 3: As a product user, I need to be able to locate the information I inputted to review my fitness and eating habits

- a. Did not finish in this sprint

F. User Stories finished ahead of time:

- a. As a product user, I want to accurately calculate my BMI so I can track any changes in my body mass.
- b. As a product user, I want to review my BMI to verify that my body mass corresponds to my glucose levels.
- **Scenario: (test done by Nooran Salim)**
 - 1. Run Sugar Analysis app from Android Studio
 - 2. Emulator opens on home page
 - 3. Click on BMI tab
 - a. Enter height (Example is 70)
 - i. Units were not implemented at this point
 - b. Enter weight (Example is 200)
 - i. Units were not implemented at this point
 - c. Click calculate to get results (Example gives us 28; overweight)

Sprint Three:

G. User Story One: As a product user, I need the notification to transfer me to the homepage of the app so I can view my blood sugar check schedule.

- a. Did not finish in this sprint

H. User Story Two: As a product user, I need to be able to locate the information I inputted to review my fitness and eating habits.

- **Scenario: (test done by Nicole Ng)**
 - a. Run Sugar Analysis app from Android Studio
 - b. Open Emulator to home
 - c. Main page: Enter blood sugar amount for each period
 - d. Check the inputted amounts in:
 - Progress Tab: calculated average blood sugar intake per day and displays in weekly bar chart
 - Log Tab: blood sugar amount per 'period' (i.e. Before Breakfast, After Lunch, etc.) listed daily in list format along with inputted weight and height
- **Scenario: (test done by Nikita Thumma)**
 - a. Run Sugar Analysis app from Android Studio
 - i. After 10 seconds it opens the Emulator to the Home Page also known as the personal information page

- b. In the Home Page, Enter your blood sugar levels (Before Dinner, After Dinner, Before Workout, After Workout, Before Breakfast, After Breakfast, Before Lunch, After Lunch) in each box
- c. Check progress in the Weekly Progress Tab (compares your blood sugar levels on a nice organized graph)
- d. The Logging Tab will list the date, time, height, and weight that corresponds to the specific period chosen in an organized list view
 - i. Ex: Before Lunch: Fri
Time: 02:41:47 PM
Weight: 200
Height: 6

Sprint Four:

I. User Story One: Users can set notifications to check blood sugar so they can constantly track their glucose.

J. Scenario: (test done by Nikita Thumma)

1. Run Sugar Analysis app from Android Studio
 - a. After 10 seconds it opens the Emulator to the Home Page (personal information page)
2. Click on the Notification Alerts Tab at the bottom
 - a. Click the white box to pull up the keyboard
 - b. Enter in Military Time you want notification to be sent
 - i. (13:43 PM)
 - ii. Exit out of the keyboard by clicking on the bottom arrow
 - c. Click corresponding label for notification (Before Dinner, After Dinner, Before Workout, After Workout, Before Breakfast, After Breakfast, Before Lunch, After Lunch)
 - d. Notification pops up at selected time (13:43PM) and says “Blood Sugar Tracker: Please Check Blood Sugar”
 - i. Redirects to the Home page once notification is clicked where you can enter your blood sugar in the corresponding label

- Scenario:

1. Run the Sugar Analysis App from Android Studio
2. Emulator will open to Home Page
3. Click on the Alerts Icon Tab
 - a. Click on the Text Box to pull up keyboard
 - b. Enter in a time that user desires to be notified

- i. This is military time.
- c. Set the time-of-day label for the entered time by clicking on the button with the label name.
- d. Wait for the notification to pop up at the entered time, with the according label.

K. User Story Two: Users can view their medical contacts and call their doctor in case they need medical attention asap. (Add-on feature)

- **Scenario: (test done by Nooran Salim)**

1. Run Sugar Analysis app from Android Studio
2. Emulator opens on home page
 - a. Enter blood sugar into textbox (Example is 180 for Before Breakfast)
 - b. Click on corresponding circle to send info to the database
3. Click on progress tab to switch view
 - a. View bar chart with progress logged for the day (Example is 180 on Sunday)
4. Click on BMI tab to switch view
 - a. Enter which units you're using, SI or US (Example is US)
 - b. Enter height (Example is 68 inches)
 - c. Enter weight (Example is 150 pounds)
 - d. Click calculate to get results (Example gives us 22; normal)
5. Click on the log tab to view record of log
 - a. Display shows the following:
 - i. Before Breakfast Blood Sugar: 180
 - ii. Date: Sun
 - iii. Time: 03:52:22 PM
 - iv. Height: 68
 - v. Weight: 150
6. Click on Notification tab
 - a. Set Military time of when you want notification (Example is 15:56)
 - b. Click on the time of day you want to set notification for (After Dinner example)
 - c. Notification pops up
 - i. Click on home to get redirected to home page
7. Click on extras tab
 - a. Read information on blood sugar level meanings
 - b. Phone number provided for Diabetes Association

8. Click on Settings
 - a. Can view/edit input from test in Sprint 1
9. Test complete; app is functioning